ADLER
ORTHO

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ADLER ORTHO

INSTRUCTIONS FOR HIP ARTICULAR PROTHESES - PANTHEON PFR

GENERAL PRODUCT INFORMATION

Even though the surgeon can now count on the means necessary to restore mobility, correct deformity, and reduce pain for many patients, thanks to the advancements in partial and total joint replacement, which generally allow him to achieve these goals successfully, it must be remembered that the prostheses used are manufactured from metal, ceramic and plastic materials and that any joint replacement system, therefore, cannot be expected to withstand activity levels and loads to avoid normal healthy bone. In addition, the system will not have the same strength, reliability or duration as a natural human joint.

- In using total joint prostheses, the surgeon should be aware of the following:
 - The correct selection of the prosthesis is extremely important.
 - The potential for success in total joint replacement is increased by selection of the proper size, shape, and design of the prosthesis. Total joint prostheses require careful seating and adequate bone support. Smaller sized implants are intended for patients with small bone and normally slight weight. Such components could be inappropriate for other patients. Surgeons are encouraged to use their best medical judgement when choosing the proper implant size regardless of the endosteal area of the bone.
 - In selecting patients for total joint replacements, the following factors can be critical to the eventual success of the procedure.
 - Patient's weight. An overweight or obese patient can produce high loads on the prosthesis, which can lead to failure of the prosthesis. This becomes a major consideration when the patient is small boned and a small size prosthesis must be used.
 - Patient's occupation or activity. If the patient is involved in an occupation or activity, which calls for substantial walking, running, lifting, or muscle strain, the resultant forces can lead to failure of the fixation, the device or both. The prosthesis will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic functional expectations.
 - Condition of spinal, mental, illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the prosthesis, leading to failure or other complications.
 - Foreign body sensitivity. In case of suspected material sensitivity, appropriate tests should be made prior to material selection or implantation.

Caution: Patients with renal insufficiency may be sensitive to potential metal ion release. Also, since little is known about the support of potential metal ion release across the placenta, these devices should be used with caution in women of childbearing age.

GENERAL

The biocompatible materials used in the manufacturing of implantable medical devices show an excellent resistance to corrosion and conform to specific international standards:

- Stainless steel, standard ISO 5832 - 1 or 9, symbol: Stainless steel
- Titanium alloy, standard ISO 5832 - 3, symbol: Ti6Al4V
- Cobalt-Chrome-Molibdenum alloy, standard ISO 5832 - 4 or 12, symbol: Co-Cr-Mo

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- Alumina ceramics, standard ISO 6474, symbol: Al2O3
- Ultra high molecular weight polyethylene, standard ISO 5834 – 1 and 2, symbol: UHMWPE
- Unalloyed titanium, standard ISO 5832 - 2, symbol: Ti. The symbol of the material used to manufacture the device is specified on the product label.

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The various components of implantable medical devices manufactured by Adler Ortho are designed in such a way as to ensure their interchangeability and compatibility.

- In particular:
 - if ceramic heads are intended to be used, modular necks and mono-block stems manufactured by Adler Ortho should only be employed in combination with ceramic heads manufactured by Adler Ortho. The company performs specific tests on these heads, to guarantee a safe coupling to the Adler Ortho taper.
 - Metal / Metal inserts have to be used only in combination with Metal / Metal heads manufactured by Adler Ortho.
- Adler Ortho take no responsibility for the association of their implants with components from another manufacturer.

INAPPROPRIATE COMBINATIONS OF MATERIALS IN ARTICULATING SURFACES

- Stainless steels / titanium alloys
- Stainless steels / stainless steels
- Stainless steels /unalloyed titanium
- Stainless steel (ISO 5832 - 1) / cobalt chrome alloys
- Unalloyed titanium / unalloyed titanium
- Unalloyed titanium / titanium alloys
- Unalloyed titanium / cobalt chrome alloys
- Unalloyed titanium / ultra high molecular weight polyethylene
- Titanium alloys / cobalt chrome alloys
- Ceramics/metal
- Zirconia ceramics /alumina ceramics

INAPPROPRIATE COMBINATIONS OF METALS IN NON ARTICULATING CONTACT SURFACES

- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / cobalt chrome alloys
- Stainless steel (excluding the stainless steel described in ISO 5832 - 9) / unalloyed titanium

INDICATIONS

- Complex Primary and Revision surgery such as
 - Deformities, such as dysplasia with complex femoral reconstruction and leg lengthening
 - Fractures and necrosis of bone with excessive bone loss in the proximal femoral region
 - Revisions, where other devices or treatments have failed
 - Primary oncology
 - Metastatic diseases

CONTRA-INDICATIONS

Contra-indications may be relative or absolute. The articular problems must be evaluated case by case, taking into account alternative surgical options (osteotomy, arthrodesis, amputation, etc.).

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The following examples are regarded as contra-indications:

- infection, septicæmia, and osteomyelitis constitute cases of absolute contra-indication;
- Pre-existing infection with or without loosening
- Movement of one or several prosthetic elements due to mechanical overloading, osteoporosis, etc.
- Dislocation of the prosthesis or fracture of the bone due to traumatism from physiotherapy
- Extra-articular pathology: pleibothrombosis, pulmonary embolism, etc.
- fatigue fracture of prosthetic components can occur as a result of : trauma, strenuous activity, improper alignment or duration of service.

Given their complexity, it is advised that, in the selection of an articular implant, the surgeon studies the clinical indications and contra-indications of that particular implant and that the surgeon uses his personal experience and judgement .

POSSIBLE COMPLICATIONS, SIDE-EFFCTS

Apart from possible per-operation complications, a prosthetic implant can also be subject to:

- Pre-existing infection with or without loosening
- Movement of one or several prosthetic elements due to mechanical overloading, osteoporosis, etc.
- Dislocation of the prosthesis or fracture of the bone due to traumatism from physiotherapy
- Extra-articular pathology: pleibothrombosis, pulmonary embolism, etc.
- fatigue fracture of prosthetic components can occur as a result of : trauma, strenuous activity, improper alignment or duration of service.

VIGILANCE SYSTEM (EUROPEAN DIRECTIVE 93/42/EEC)

The users are party to the setting up of the vigilance system after introduction of medical device on the market.

They are obliged by law to report to the competent authorities as well as the manufacturer any incident concerning the medical device liable to cause or having caused serious risk for the patient or user health.

VERIFICATION PRIOR TO USE - PRECAUTIONS

It is essential to have all the sizes listed in the catalogue in the operating theatre prior to the surgery to ensure the optimum size for the patient is selected.

Presentation cases and instruments are available to achieve this. Maintenance and handling must be carried out in accordance with the recommendations of standard ISO 8828.

Handling must be limited and performed with extreme care (non-powered gloves must be worn for the manipulation of devices coated with hydroxyapatite) so as to avoid damage/scratching, shock or contact with materials liable to alter the surface of the implant). In all cases the final confirmation of the condition of the implant lies with the surgeon. Special care must be taken with ceramic devices which are not to be re-used even when in the absence of any apparent damage. Only use a polyacetal tip to introduce ceramic devices.

All implantation all foreign bodies such as pieces of tissue and bone or cement particles must be removed. They could cause abnormal wear of the articulating surfaces.

If more layers are used to assemble the implants, confirmation of the compatibility of the layer angles must be made. Before the neck is inserted, the female cone of the stem should be thoroughly cleaned. All components should be thoroughly

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cleaned before their assembly.

Use the femoral prostheses for dysplasia only in special cases, taking into account that the body weight may affect the mechanical resistance of the prosthesis considerably. Under no circumstances should a damaged or scratched prosthesis or a prosthesis having undergone any inappropriate or unauthorised treatments be implanted, this includes implants that have already been used. Never steam sterilize/sterilize: femoral prostheses with ceramic femoral heads sealed on the stem, ceramic, plastic and/or metal/plastic components. Please check the sterilization by verifying the colour of the sterilization indicator (red in case of gamma sterilization, green in case of ethylene oxide sterilization). Please verify the vacuum tightness of the packaging of the metallic devices. Do not use a device in which the packaging has been damaged or is questionable because it could call into question the integrity of the packaging and the sterility of the device. Please verify the sterilization expiration date.

RE-USE

Implants must never be reused as their mechanical and/or functional characteristics could be definitively compromised so their service life could be reduced.

INSTRUMENTS SUPPLIED WITH THE PROSTHESIS

Instruments supplied with the prosthesis must be sterilized by means of a validated and controlled sterilization method.

For metal instruments we recommend steam sterilization (Autoclave) usually employed in hospitals.

134°C cycles for 18 minutes do not cause any deformation or damage to the plastic or metal instruments.

It is important to avoid any contact with objects or corrosive substances (e.g.: strong acid solutions, bases, solvents) which could adversely alter the material.

DEVICES SUPPLIED IN A STERILE STATE

The metallic devices supplied by Adler Ortho in a sterile state are sterilized in vacuum by means of gamma rays radiation at a nominal dose of 25 kGy. The polyethylene devices supplied by Adler Ortho in a sterile state are sterilized by means of ethylene oxide.

The label clearly shows the sterilisation indicator and the indication "STERILE".

RESTERILISATION

All implants supplied by Adler Ortho must not be re-sterilised by the purchaser. Adler Ortho does not take any responsibilities for the use of implants re-sterilised by the purchaser.

STORAGE CONDITIONS

The implants must be stored under dry conditions and be protected from sunlight, extremes in temperature and contamination by particles.

PATIENT INFORMATION

The patient destined to receive the implant must be informed that the security and durability of the implant depends also on his/her behaviour, particularly as concerning his/her activity and body weight. For any further information please contact your supplier.

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